


### **Quick Start Manual PLENTY Vaporizer**

## STORZ & BICKEL

1. Open Filling Chamber




**2.** Grind plant material




3. Fill Filling Chamber




Always fill the Filling Chamber completely up to the upper edge







#### correct




Completely filled Filling Chamber






Filling Chamber only partially filled




**4.** Switch on the device (orange power switch) and press auto off release handle.



**5.** Vaporization temperature scale



**6.** Choose temperature



Refer to the Instructions for Use for important safety and user advice. Do not operate the device without first reading and understanding the Instructions for Use!


# © by STORZ & BICKEL GmbH · PAL-30-092 12-2021 · Subject to alterations · All rights reserved

## plenty

### Disassembling, Cleaning, Assembling




**1.** Expanded view of the Vaporization Unit



Free Shipping Worldwide



**2.** Pull the Mouthpiece away from the Cooling Coil



**3.** Pull the Cooling Coil out of the Filling Chamber Cap




**4.** Push out the Upper Screen with help of the Cleaning Brush



5. Disassembling the Filling Chamber



**6.** Push out the Filling Chamber Cylinder



7. Push out the Lower Screen



**8.** Manual cleaning with alcohol (min. 90%) or dishwashing soap



**9.** Filling Chamber Cylinder, Filling Chamber Housing, positioning surfaces

